

 Read this entire package insert carefully prior to use.

 Federal law (U.S.A.) restricts this allograft bio-implant for use by a licensed clinician only.

DESCRIPTION

This allograft bio-implant was processed from donated human tissue, resulting from the generous gift of an individual or his/her family. The bio-implant was processed per the surgeon's specifications, packaged in a nutrient medium with added antibiotics, Gentamicin and Vancomycin, and then refrigerated. Processing was performed under aseptic conditions and is culture negative.

INDICATIONS FOR USE

This allograft bio-implant is intended for the repair or replacement of damaged osteochondral tissue.

CONTRAINDICATIONS

The contraindications include, but are not limited to:

- Use in any patient who has a known or suspected allergy to any of the antibiotics and/or processing reagents listed in this package insert.

WARNINGS AND PRECAUTIONS

The same medical/surgical conditions or complications that apply to any surgical procedure may occur during or following implantation. The surgeon is responsible for informing the patient of the risks associated with their treatment and the possibility of complications or adverse reactions. As with any allograft bio-implant, the potential for transmission of infectious agents exists.

This bio-implant may contain residuals of antibiotics (Gentamicin and Vancomycin) and human proteins (insulin, transferrin, human serum albumin). Caution should be exercised if the patient has a known sensitivity to any of these residual elements.

STORAGE REQUIREMENTS

The distributor, intermediary and/or end-user clinician or facility is responsible for storing this allograft bio-implant under appropriate conditions prior to further distribution or implantation. Bio-implants must be stored as listed in the table below.

DO NOT FREEZE.

Temporary excursions above freezing (0°C to 1°C) may occur during shipment.

Preservation Method	Storage Temperature
Refrigerated and stored in a nutrient medium	Store between 1°C and 10°C.

POTENTIAL ADVERSE EVENTS

Potential adverse events or outcomes include, but are not limited to, disease transmission, infection, allograft tissue rejection, allergic reaction to residual processing reagents, re-operation and/or death.

Promptly report any adverse event(s) or outcome(s) potentially attributable to the allograft bio-implant (See **COMPLAINTS AND RETURNS** section).

DONOR SCREENING AND TESTING

All donors have been screened and tissues recovered, processed, stored, tested and distributed in accordance with current U.S. federal regulations as promulgated in 21 CFR 1270 and 1271, current Standards for Tissue Banking set forth by the American Association of Tissue Banks (AATB) and international laws and regulations as required.

This allograft bio-implant was deemed suitable for implantation by LifeNet Health. A physician medical director evaluated the following donor variables to determine donor suitability: infectious disease test results, current donor medical history, behavioral risk assessment interview, physical assessment, relevant medical records, including previous medical history, laboratory test results, and autopsy or coroner reports (if performed).

All donors are tested for relevant infectious diseases. Testing is performed by laboratories that are registered with the U.S. Food and Drug Administration (FDA) and certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR 493. Test methods that are FDA-licensed, approved, or cleared for donor screening are used as available. The following test criteria were met for the donor of this allograft bio-implant:

Required Infectious Disease Testing	
Test	Acceptance Criteria
HBcAb: Hepatitis B Total Core Antibody	Negative/Non-Reactive
HBsAg: Hepatitis B Surface Antigen	Negative/Non-Reactive
HCV NAT: Hepatitis C Virus Nucleic Acid Test	Negative/Non-Reactive
HCVAb: Hepatitis C Antibody	Negative/Non-Reactive
HBV NAT: Hepatitis B Virus Nucleic Acid Test**	Negative/Non-Reactive
HIV-1 NAT: Human Immunodeficiency Virus Type 1 Nucleic Acid Test	Negative/Non-Reactive
HIV 1/2 Ab: Human Immunodeficiency Virus Types 1/2 Antibody	Negative/Non-Reactive
RPR/STS or Equivalent: Syphilis	Confirmatory Negative/Non-Reactive
HTLV I/II Ab: Human T-Lymphotropic Virus Types I/II Antibody*	Negative/Non-Reactive

*Not required for donors recovered after March 31, 2010. Performed as required by international laws and regulations.

**Not required for donors recovered prior to December 16, 2016. Performed as required by international laws and regulations.

SEE BACK FOR GRAFT PREPARATION INSTRUCTIONS

INSTRUCTIONS FOR USE



It is important to read and understand the following instructions prior to clinical use. Improper preparation technique may adversely affect handling properties and/or performance.

GENERAL INSTRUCTIONS:

- Use on a single occasion for a single patient only.
- Once the packaging is opened, the bio-implant must be used for the current procedure or discarded.
- Inspect the bio-implant, inner and outer packaging, and labels carefully:
 - Do not use past the expiration date as indicated on the label.
 - Do not use if the bio-implant is damaged or the packaging integrity is compromised.
 - If liquid is present inside the outer tray, examine the inner tray seal for integrity. Do not use the bio-implant if the inner tray seal is found to be compromised.
 - Do not use if there are discrepancies in label information.
- Use aseptic technique at all times.
- Do not sterilize.
- Keep the bio-implant stored according to recommended storage instructions until preparing it for implantation.
- The bio-implant should be handled gently to protect the cartilage.

PREPARATIONS FOR USE:

Circulator:

- This graft has been processed aseptically and packaged in a dual tray sterile packaging system.
- Remove the dual tray from the product box (if present).
- Using sterile technique, peel open the foil lid on the outer tray by pulling on the peel tab located on a corner of the tray. Do not touch the inside surface of the outer tray, or the inner tray itself while opening the packaging. The inside of the outer tray and inner tray is sterile.
- Firmly hold the outer tray below the top surface, and present the inner tray to the sterile scrub technician for removal.
- During transfer, the sterile technician should not touch the outside of the outer tray packaging as this is not sterile.

Sterile Member:

- Being careful not to touch the outside of the outer tray, remove the inner tray containing the bio-implant.
- Place the inner tray on a surface in the sterile field. Open the tray by peeling back the top foil surface by using the peel tab located on a corner of the tray.
- The Fresh graft may be kept in the opened inner tray until ready for use in surgery as the transport medium is optimized for cartilage preservation.
- Fresh bio-implants should be handled gently and kept moist at all times, either soaking in sterile media or covered with wet gauze.
- Lavage the Fresh bio-implant prior to implantation.

TRACEABILITY

It is the responsibility of the end-user to maintain recipient records for the purpose of tracking tissue post-implantation. As a courtesy to the end-user clinician or facility, LifeNet Health has enclosed a Graft Implant Tracking Card to assist in the post-implantation tracking. Please refer to the enclosed card for additional instructions.

COMPLAINTS AND RETURNS

For further information on returns or to report a complaint or adverse event, please contact your authorized distributor or LifeNet Health Client Services (available 24 hours a day) at 1-888-847-7831 (inside the U.S.) or 00+1-757-464-4761 ext. 2000 (outside the U.S.) and have the bio-implant's identification number available (see label).

WARRANTY STATEMENT

Due to the inherent variability of allograft tissue, biological and biomechanical properties cannot be guaranteed by LifeNet Health.

Manufactured by:

LifeNet Health
1864 Concert Drive
Virginia Beach, VA 23453

Source Establishment: LifeNet Health CTO #100038

63-0238.01

Fresh allografts are covered by one or more of the following US patents:
USD793,251

The LifeNet Health logo is a registered trademark of LifeNet Health.
©2018 LifeNet Health. All rights reserved.

North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ Asia □ Middle East
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu