



Preparing the Autograft for Shipment

Using the Recovery Kit

OSTEOCLEANSE

A turn-key solution for the management of autologous tissue.

OUTSIDE THE STERILE FIELD

1



1. Open the box by cutting the tape on the top and sides of the box.

2



2. Remove the Ziploc® bag containing all the forms.

3



3. Remove the Sterile Kit.

4



4. Using sterile technique unwrap the Kimguard® wrap and pass the inner contents on to the sterile field.

5



5. (Sterile Person) Take the kit and remove the components from within the Tupperware® container.

6



6. After the craniectomy, remove any hardware attached to the explanted bone flap.

7



7. Debride the bone flap of soft tissue.

8



8. Rinse with 2–3 liters of sterile solution.

9



9. Place the bone flap inside the 9"x12" Ziploc® bag.

{continued inside}

Preparing the Autograft for Shipment (continued)

SAFE & BIOCOMPATIBLE


CLINICALLY EFFECTIVE

COST EFFECTIVE



10

10. Fold the Ziploc® bag over several times and seal.



11

11. Wrap the Ziploc® containing the bone flap with the Kimguard® wrap (use only one of the two provided).



12

12. Fold and place the second Kimguard® inside the Tupperware® container.



13

13. Place the bone flap inside and securely fasten the lid on the Tupperware® container.



14

14. Use the included Sharpie® pen to write:

- Patient's name
- Medical record number
- Date of explanation
- Time of explanation



15

15. Place the sealed Tupperware® container inside the provided sterile plastic bag (1 of 2) and knot the bag or close it with a plastic tie.



16

16. Place the bagged Tupperware® container inside the provided sterile plastic bag (2 of 2) and knot or close it with a plastic tie.



17

17. Pass the double bagged Tupperware® container off the sterile field.

Preparing the Kit for Shipment



1. Complete the information on the white tag (can be found inside the literature pack in the kit) and attach it to the bagged Tupperware® container with the bone flap using a locking plastic tie.

a. Do not fill in the Donor ID field.



2. Spread the plastic liner bag in the insulated shipper and fill the box halfway with ice. Place the bagged Tupperware® container (lid side up) on top of the ice.

Note: LifeNet Health cannot accept autologous tissue shipments if the information provided is incomplete. This information is required by the American Association of Tissue Bank Standards for Tissue Banking.




3. Fill the remainder of the insulated shipper with ice, making sure that the bagged container holding the autograft is completely surrounded by ice and the insulated shipper is filled (leave about 1" for the lid).



4. Knot the plastic liner bag.



5. Place the lid on the insulated shipper.



6. Complete included Request Form; fax to LNH at 1-800-617-1389 or email to osteo cleanse@lifenethealth.org.



7. Place a copy of the completed OsteoCleanse Service Request Form in the re-sealable (Ziploc®) plastic bag. Place the bag on top of the container lid.



8. Seal the cardboard box with packing tape. If the flap is to be cleaned and disinfected continue to step 9, if the flap is for storage only, please call Client Services at 1-888-689-5749 for shipping instructions.



9. Write your hospital's address in the appropriate field on the provided airway bill and place it in the clear sleeve. **Do not seal the sleeve.**

Flaps for storage only should be shipped to: LifeNet Health, 1864 Concert Drive Virginia Beach, VA 23453.

{continued on back}

Preparing the Kit for Shipment (continued)



10

10. Attach the sleeve with the airway bill to the top of the box. **Do not seal the sleeve.**



11

11. On the fluorescent pink shipping label located on the side of the cardboard box, write the:

- Name of the person who packaged the shipment
- Medical Record Number
- Current date
- Current time



12

12. Call LifeNet Health's OsteoCleanse line at 1-888-689-5749 to arrange for pick-up.

NEXT STEPS: WHAT TO EXPECT

1



1. Upon the graft's arrival at LifeNet Health, you will receive an e-mail confirming receipt of the graft.

2



2. Once the OsteoCleanse process is complete you will receive an e-mail notifying you that the process is complete and ready for shipment back to your facility, or it can be stored at LifeNet Health.

Contact

Toll Free: (888) 689-5749

Fax: (800) 617-1389

Email: osteocleanse@lifenethealth.org

24
hours
per day

7
days
per week

365
days
per year

WHEN THE GRAFT ARRIVES BACK AT YOUR FACILITY

1



1. Inside the shipping container you will find the graft. See thawing directions and attached Instructions for Use (IFU).

2



2. If the graft arrives in a reusable shipper, inside you will find a pre-paid return label. Place label on box for shipment back to LifeNet Health.

Room Temperature (Standing Air)	Warm Solution Soak (37°C–42°C)	Room Temperature Solution Soak
20 minutes	3 minutes	5 minutes

LifeNet Health
1864 Concert Drive
Virginia Beach, VA 23453

www.LifeNetHealth.org

