



OraGRAFT[®]

Cancellous Blocks and Cubes

Clinical Overview The OraGraft product line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects.

- Applications**
- Ridge Maintenance (Horizontal & Vertical Augmentation)
 - Craniofacial Reconstruction
-

- Features & Benefits**
- **Increased Safety:** Sterile with an SAL 10⁻⁶.
 - **Convenient Storage:** Ambient temperature with five-year shelf life.
 - **Ease of Use:** Can be affixed to surgical site using standard bone screws
 - **Rapid Remodeling:** Solid cancellous bone allows for space maintenance and rapid remodeling due to trabecular architecture.
 - **Match to Defect:** Easy to trim, can be shaped using standard instrumentation.
-





OraGraft Cancellous Blocks and Cubes

Ambient Storage*/5 Year Shelf Life

Size	Order Code
10 x 10 x 10 mm	CANCUBE1
15 x 15 x 15 mm	CCUBE-01
15 x 15 x 8 mm	CCUBE-02
15 x 30 x 8 mm	CCUBE-03

*While ambient room temperature has not been defined by regulatory bodies, LifeNet Health would recommend storage at 2°C to 37°C with excursions of less than 24 hours up to 40°C. If an excursion outside this range occurs, please contact LifeNet Health.

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)

