



OraGRAFT[®]

Bi-Cortical Ilium Strips

Clinical Overview The OraGraft product line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects.

- Applications**
- Ridge Maintenance
 - Craniofacial Reconstruction
-

- Features & Benefits**
- **Increased Safety:** Sterile with an SAL 10⁻⁶.
 - **Convenient Storage:** Ambient temperature with five-year shelf life.
 - **Easy of Use:** Can be affixed to surgical site using standard bone screws, easily trimmed using standard instrumentation to expose cancellous located within cortical plates.
-





OraGraft Bi-Cortical Ilium Strips

Ambient Temperature*/5 Year Shelf Life

Size	Freeze-Dried	Preservon
15 x 20 mm	IS SML	ISP SML
15 x 30 mm	IS MED	ISP MED
20 x 40 mm	IS	ISP

*While ambient room temperature has not been defined by regulatory bodies, LifeNet Health would recommend storage at 2°C to 37°C with excursions of less than 24 hours up to 40°C. If an excursion outside this range occurs, please contact LifeNet Health.

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)

