



FlexiGRAFT[®]


Quadriciceps Tendons

Clinical Overview FlexiGraft Quadriciceps Tendons can be used for ACL and PCL reconstruction procedures.

- Applications**
- ACL Reconstruction
 - PCL Reconstruction

- Features & Benefits**
- **Convenience:** No time-consuming harvest required. Pre-shaped grafts eliminate prep time and further increase operating room efficiency.
 - **Patient-Friendly:** Eliminates donor site morbidity and associated pain from autograft recovery. This makes the procedure less invasive and potentially decreases OR time. Less OR time can mean less time under anesthesia and less tourniquet time.¹
 - **Sterile:** Achilles tendons are sterilized using proprietary and patented Allowash XG[®] technology. This provides a sterility assurance level (SAL) of 10^{-6} , without compromising the construct's inherent biomechanical properties.²





FlexiGraft Quadriceps Tendons

Frozen Storage (-40°C to -80°C) / 5 Year Shelf Life

Description	Tendon Length (TL)	Bone Block Size	Order Code
Quadriceps Tendon with Bone Block	65 - 95 mm	Width = 10 - 12 mm (W) Length = \geq 20 mm (L)	FQUADB

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)

References

- Oro et al. Autograft Versus Allograft: An Economic Cost Comparison of Anterior Cruciate Ligament Reconstruction. *Arthroscopy* 2011; 27(9):1219-1225.
- LifeNet Health ES-04-015.

