



OraGRAFT[®]

Ribs

Clinical Overview The OraGraft allograft line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects

- Applications**
- Ridge Maintenance
 - Craniofacial Reconstruction

- Features & Benefits**
- **Highest Safety Level:** Sterile with an SAL of 10^{-6} .
 - **Convenient Storage:** Ambient temperature with a 5 year shelf life.
 - **Ease of Use:** Can be affixed to surgical site using standard bone screws.
 - **Versatile:** Solid cortical shell with cancellous contained within.
 - **Match to Defect:** Easy to trim, can be shaped using standard instrumentation.



OraGraft Ribs

Ambient Temperature*/5 year shelf life

Sizes	Order Code
60 - 115 mm	RIB

*While ambient room temperature has not been defined by regulatory bodies, LifeNet Health would recommend storage at 2°C to 37°C with excursions of less than 24 hours up to 40°C. If an excursion outside this range occurs, please contact LifeNet Health.

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)