

## FREQUENTLY ASKED QUESTIONS

### Decellularized Dermis Packaging Summary

LifeNet Health's line of decellularized dermis grafts (OrACELL<sup>®</sup>, DermACELL<sup>®</sup> and ArthroFlex<sup>®</sup>) are packaged in double sterile barrier Aluminum Oxide Polyester Pouches. These pouches are packaged within an outer card board sleeve. These products are shipped in insulated shippers with cooling bricks which have been validated to maintain adequate temperatures for 72 hours of distribution using the ISTA (International Safe Transit Association) 7D winter and a modified ISTA 7D summer profile. The temperature profiles tested were more robust than the 7D standard and address potential worst case scenarios for summer shipments. The packaging contains a temperature sensitive dot that will turn from white to pink or red if the upper temperature limit has been exceeded. The decellularized dermis should not be used if the temperature dot appears to be a color other than white.