



ArthroFLEX[®]

Acellular Dermal Matrix

Clinical Overview

- ArthroFlex is an acellular dermal matrix used to supplement soft tissue repairs such as the Achilles tendon and rotator cuff. It has $\geq 97\%$ DNA removal, comes ready-to-use right out of the package, can be stored at room temperature (15°C - 30°C), and has a Sterility Assurance Level (SAL) of 10^{-6} , which is equal to medical device-grade sterility

Applications

- Augment repairs of the rotator cuff
- Augment repairs of the Achilles, biceps, quad and patellar tendons
- Reinforce peroneal/tibialis tendons and flexor/extensor tendons of hand and foot
- Reinforce capsular repairs
- Suture reinforcement of soft tissue repairs

Why Use

- **Improved Clinical Outcomes:** ArthroFlex demonstrated decreased re-tear rates and improved patient reported outcome scores after one and two years^{2,3}
- **Increased Safety:** Sterile with a SAL of 10^{-6} means lower risk of disease transmission
- **Decellularized:** Uses Matracell[®] technology to remove $\geq 97\%$ DNA and cellular remnants to decrease likelihood of immune response^{1,4}
- **Promotes Rapid Healing:** Biocompatible and retains natural growth factors, collagen and elastin^{1,4}
- **Convenient:** Preservon[®] technology allows graft to be stored fully hydrated at room temperature (15°C - 30°C) so it can be used right out of the package, without lengthy OR reconstitution processes^{1,4}
- **Versatility:** Available in multiple sizes and thicknesses





ARTHROFLEX ACELLULAR DERMAL MATRIX

Room Temperature Storage (15°C - 30°C) / 3 Year Shelf Life

Size	Thickness	Order Code
35 x 35 mm	1.5 mm (1.26 - 1.74)	AFLEX100
40 x 70 mm	1.5 mm (1.26 - 1.74)	AFLEX101
50 x 90 mm	1.5 mm (1.26 - 1.74)	AFLEX103
15 x 140 mm	1.5 mm (1.26 - 1.74)	AFLEX150
35 x 35 mm	2.0 mm (1.76 - 2.25)	AFLEX200
40 x 70 mm	2.0 mm (1.76 - 2.25)	AFLEX201
25 x 30 mm	2.0 mm (1.76-2.25)	AFLEX 202
40 x 50 mm	3.0 mm (2.5 - 3.5)	AFLEX300
40 x 70 mm	3.0 mm (2.5 - 3.5)	AFLEX301
50 x 60 mm	3.0 mm (2.5 - 3.5)	AFLEX302
40 x 40 mm	1.0 mm (1.26 - 1.74)	AFLEX400
40 x 70 mm	1.0 mm (1.26 - 1.74)	AFLEX401
30 x 40 mm	0.5 mm (0.3 - 1.0)	AFLEX500
BioWasher® 10 x 14 mm, 2 pack	2.0 mm (1.76 - 2.24)	AFLEX822

References

1. Data on file at LifeNet Health TR-0272 to 65-0272
2. Gilot GJ, Alvarez AM, Barcksdale L. Outcome of large to massive rotator cuff tears repaired with and without extracellular matrix augmentation: A prospective comparative study. *Arthroscopy*. 2015 Apr 17. pii: S0749-8063(15)00141-3.
3. Denard PJ, Brady PC, Adams CR, Tokish JM, Burkhart SS. Preliminary Results of Arthroscopic Superior Capsule Reconstruction with Dermal Allograft. *Arthroscopy*. 2018;34(1):93-9.
4. Moore MA, Samsell B, Wallis G, Triplett, S, Chen S, Linthorst Jones A, Quiz X. "Decellularization of Human Dermis Using Non-Denaturing Anionic Detergent and Endonuclease: A Review." *Journal of Cell and Tissue Banking*.

68-60-081-02.02

LifeNet Health, the LifeNet Health logo, ArthroFlex, BioWasher, Preservon and Matracell are registered trademarks of LifeNet Health. The Arthrex logo is a registered trademark of Arthrex, Inc. ©2020 LifeNet Health. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu