



FlexiGRAFT[®]

Meniscus


Applications

- Meniscal transplant
-

Why Use

- **Convenience:** Processed with tibial bone block to allow for surgical technique flexibility such as double bone plug technique, keyhole technique or dove tail technique.
- **Safety:** Menisci are processed using proprietary and patented Allowash[®] technology. Passes USP Sterility test prior to release
- **Quality:** All of LifeNet Health's menisci are processed from donors 45 years old or younger to ensure best possible meniscus quality
- **Custom Sizing:** Our extensive measurement parameters are designed to find the most precise fit for the patient
- **Custom Matching:** Upon receipt of appropriate films, x-rays or MRI and using the "Pollard" sizing technique, our technical staff will provide the surgeons with precise measurements and graft options





FlexiGraft Meniscus

Frozen storage / 5 year shelf life

Description	Minimum Bone Height & Width Others "Custom"	Order Code
Right Medial Meniscus	15 mm (in 1 mm increment)	FMN RM
Right Lateral Meniscus	15 mm (in 1 mm increment)	FMN RL
Left Medial Meniscus	15 mm (in 1 mm increment)	FMN LM
Left Lateral Meniscus	15 mm (in 1 mm increment)	FMN LL

Measurements taken and provided on dissection sheet and offered to surgeon:

Meniscus

- (A)** Length (27 to 55 mm)
- (B)** Width (25 to 45 mm)
- (C)** Height (\geq 3 mm)
- (D)** Inferior Rim Width

Tibial Bone Block

- (E)** Length
- (F)** Width (\geq 15 mm)
- (G)** Height (\geq 15 mm)

Tibial bone block width minimum = 15 mm (in 1 mm increments)

Tibial bone block height minimum = 15 mm (in 1 mm increments)

68-60-188.00

The LifeNet Health logo is a registered trademark of LifeNet Health.
FlexiGraft and Allowash are registered trademarks of LifeNet Health.
©2018 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu