



OraGRAFT[®]

Cancellous Mineralized Powder

Clinical Overview The OraGraft line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects.

- Applications**
- Periodontal Defects
 - Ridge Maintenance
 - Extraction Site Preservation
 - Sinus Grafting
 - Craniofacial Reconstruction
-

- Features & Benefits**
- **Increased Safety:** Sterile with an SAL 10⁻⁶.
 - **Convenient Storage:** Ambient temperature with three-year shelf life.
 - **Graft Incorporation:** Osteoconductive potential with a combination of open trabecular architecture and dense particles allows for healing through “creeping substitution.”
-





OraGraft Cancellous Mineralized Powder

Ambient Temperature*/3 year shelf life

250 - 1000 microns

Size	Order Code
0.50 cc	O CAN-0.5A
1.00 cc	O CAN-1.0A
2.00 cc	O CAN-2.0A

1000 - 2000 microns

Size	Order Code
0.50 cc	O CAN-0.5B
1.00 cc	O CAN-1.0B
2.00 cc	O CAN-2.0B

*While ambient room temperature has not been defined by regulatory bodies, LifeNet Health would recommend storage at 2°C to 37°C with excursions of less than 24 hours up to 40°C. If an excursion outside this range occurs, please contact LifeNet Health.

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)

