

ReadiGRAFT® BLX

Fibers



ReadiGRAFT BLX Fibers are a proprietary mix of demineralized cortical fibers and mineralized cancellous chips specifically engineered to maintain the osteoinductive potential and osteoconductive properties of the graft while providing surgeons with the operative flexibility to customize hydrating options to fit patient's needs.

Brand

- ReadiGraft BLX Fibers

Clinical Applications

- Any surgical application that requires bone void filler

Features And Benefits

- **Osteoconductive:** The natural osteoconductive porous architecture of cancellous bone and the increased surface area of demineralized cortical fibers provides a scaffold for cellular attachment and proliferation.
- **Osteoinductive Potential:** The cortical fibers are demineralized using LifeNet Health's patented and proprietary PAD technology. This technology carefully exposes natural growth factors trapped within cortical bone while maintaining its inherent osteoinductive potential.
- **No Carrier:** Contains no fillers or carriers commonly used in bone putties and gels.
- **Customizable Hydration:** Allows surgeons to customize hydration options to fit patient's needs by the following hydration options:
 - Platelet Rich Plasma (PRP)
 - Bone Marrow Aspirate (BMA)
 - Whole Blood
 - Saline
 - Antibiotic Solution
- **Biocompatible:** ReadiGraft BLX Fibers are comprised of 100% natural bone matrix which is remodeled naturally during the bone healing process.
- **Sterile:** Sterilized using proprietary and patented Allowash XG technology. This technology provides a sterility assurance level of 10^{-6} , without compromising the implant's inherent osteoconductivity or osteoinductive potential.



ReadiGRAFT[®] BLX

Fibers



BLX FIBERS

Volume

2.5 cc

Order No. DF-1002
Freeze Dried

5.0 cc

Order No. DF-1003
Freeze Dried

15 cc

Order No. DF-1004
Freeze Dried

30 cc

Order No. DF-1005
Freeze Dried

1864 Concert Drive
Virginia Beach, VA 23453
757-464-4761 x2000
Fax: 757-464-6439

www.AccessLifeNetHealth.org

The LifeNet Health and Readigrft logos are registered trademarks of LifeNet Health.
©2012 LifeNet Health, Virginia Beach, VA. All rights reserved.
68-40-095-01.00



U.S. Use Only

Phone: U.S. 1-888-847-7831 • 757-464-4761 x2000 • Fax: U.S. 1-888-847-7832 or 757-464-6439
www.AccessLifeNetHealth.org