



# ReadiGRAFT®

## Cortical/Cancellous Chips

### Clinical Overview

A natural mix of cortical and cancellous bone chips designed to encourage bi-phasic remodeling.

### Applications

- Any surgical application that requires bone void filler.

### Why Use

- **Osteoconductive:** Natural bone matrix facilitates cell attachment and proliferation, and vascular in-growth
- **Bi-Phasic:** Natural mix of cortical and cancellous bone facilitates bi-phasic remodeling
- **Pre-hydrated:** Allograft bio-implants featuring Preservon® are stored in a fully-hydrated state at ambient temperatures. Preservon eliminates thawing and re-hydration time, does not require freezer storage, and does not compromise the graft's inherent osteoconductive properties<sup>1</sup>
- **Sterile:** Sterilized using proprietary and patented Allowash XG® technology which provides a sterility assurance level of 10<sup>-6</sup>, without compromising the graft's inherent osteoconductive properties<sup>2</sup>
- **Versatile:** Available in grind sizes of 1-4 mm or 1-8 mm in multiple volumes to meet surgical needs
- **Absorbent:** Absorbs and retains bioactive fluids like blood, platelet rich plasma (PRP), and bone marrow aspirate (BMA)

### Alternative Grafts:

- Readigraft Cancellous Chips
- MatriGraft Femoral Head (Grinder)

### References

1. Independent sources include the Virginia Commonwealth University Medical Center and the American Association of Mechanical Engineers. Data of file at LifeNet Health.
2. Eisenlohr LM. "Allograft Tissue Sterilization Using Allowash XG®". 2007 Bio-Implants Brief.





## ReadiGraft Cortical/Cancellous Chips

Grind Size	Volume	Freeze-Dried	Preservon
1 - 8 mm	10	CC10	PCC10
	15		PCC15
	20	CC1/4	PCC1/4
	30	CC30	PCC30
	40	CC1/2	PCC1/2
	60	CC60	PCC60
	90	CC90	PCC90

  

Grind Size	Volume	Preservon
1 - 4 mm	15	PCC15 14
	30	PCC30 14
	60	PCC60 14
	90	PCC90 14

LifeNet Health □ 1864 Concert Drive □ Virginia Beach, VA 23453  
68-60-019-01.02

The LifeNet Health logo is a registered trademark of LifeNet Health, Virginia Beach, VA.  
Preservon and Allowash XG are registered trademarks of LifeNet Health.  
©2017 LifeNet Health. All rights reserved.



**North America**  
1.888.847.7831  
orders@lifenethealth.org

**Europe**  
+ 43 1 375002710  
eu\_orders@lifenethealth.eu

**Latin America □ Asia □ Middle East**  
1.757.464.4761 ext. 2000  
internat.orders@lifenethealth.org

**www.LifeNetHealth.org**  
**www.LifeNetHealth.eu**