



ReadiGRAFT®

Cancellous Cubes and Blocks

Clinical Overview

Sterile cancellous cubes and blocks designed for use in well-defined bone voids.

Applications

- Any surgical application that requires bone void filler.

Why Use

- **Osteoconductive:** Natural bone matrix facilitates cell attachment and proliferation, and vascular in-growth
- **Pre-hydrated:** Allograft bio-implants featuring Preservon® are stored in a fully-hydrated state at ambient temperatures. Preservon eliminates thawing and re-hydration time, does not require freezer storage, and does not compromise the graft's inherent osteoconductive properties¹
- **Sterile:** Sterilized using proprietary and patented Allowash XG® technology which provides a sterility assurance level of 10^{-6} , without compromising the graft's inherent osteoconductive properties²
- **Versatile:** Available in multiple sizes and volumes to meet surgical needs
- **Absorbent:** Absorbs and retains bioactive fluids like blood, platelet rich plasma (PRP), and bone marrow aspirate (BMA)

Alternative Grafts:

- Readigraft Cancellous Chips
- Readigraft Cortical/Cancellous Chips
- MatriGraft Femoral Head (Grinder)

References

1. Independent sources include the Virginia Commonwealth University Medical Center and the American Association of Mechanical Engineers. Data of file at LifeNet Health.
2. Eisenlohr LM. "Allograft Tissue Sterilization Using Allowash XG®". 2007 Bio-Implants Brief.





ReadiGraft Cancellous Cubes

Size	Volume	Freeze-Dried	Preservon
0.5 x 0.5 x 0.5 cm	5cc	CANCUBE05 S	CANCUBE05 SP
	15cc	CANCUBE15 S	CANCUBE15 SP
Size	Volume	Freeze-Dried	Preservon
1.0 x 1.0 x 1.0 cm	8 Cubes	CANCUBE1/4	PCANCUBE1/4
	12 Cubes	CANCUBE30	PCANCUBE30
	16 Cubes	CANCUBE1/2	PCANCUBE1/2



ReadiGraft Cancellous Blocks

Size	Volume	Freeze-Dried	Preservon
20 x 20 x 30 mm	1 Block	CANBLOCK	PCANBLOCK

LifeNet Health □ 1864 Concert Drive □ Virginia Beach, VA 23453
68-60-020-01 .01

The LifeNet Health logo is a registered trademark of LifeNet Health, Virginia Beach, VA.

Preservon and Allowash XG are registered trademarks of LifeNet Health.

©2017 LifeNet Health. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ Asia □ Middle East
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu