

FlexiGRAFT®

Demineralized Fibers



FlexiGRAFT® Demineralized Fibers are a new form of demineralized bone which offers excellent handling properties without the need of a carrier. FlexiGraft Demineralized Fibers are engineered to maintain their osteoinductive potential and osteoconductive properties, while providing surgeons the operative flexibility to customize hydrating options to fit patient's needs.

Brand	<ul style="list-style-type: none"> FlexiGraft Demineralized Fibers
Clinical Applications	<ul style="list-style-type: none"> Any surgical application that requires bone void filler
Features And Benefits	<ul style="list-style-type: none"> Osteoconductive: The increased surface area of demineralized cortical fibers provides an optimal scaffold for cellular attachment and proliferation. Osteoinductive Potential: The cortical fibers are demineralized using LifeNet Health's patented and proprietary PAD technology. This technology carefully exposes natural growth factors trapped within cortical bone while maintaining its inherent osteoinductive potential. No Carrier: Contains no fillers or carriers commonly used in bone putties and gels. Customizable Hydration: Allows surgeons to customize hydration options to fit patient's needs by the following hydration options: <ul style="list-style-type: none"> Platelet Rich Plasma (PRP) Bone Marrow Aspirate (BMA) Whole Blood Saline Antibiotic Solution Biocompatible: FlexiGraft Demineralized Fibers are comprised of 100% natural bone matrix which is remodeled naturally during the bone healing process. Sterile: Sterilized using proprietary and patented Allowash XG technology. This technology provides a sterility assurance level of 10^{-6}, without compromising the grafts inherent osteoconductivity or osteoinductive potential.
Alternative Products	<ul style="list-style-type: none"> FlexiGraft Demineralized Sponge Cancellous Cortical Cancellous



FlexiGRAFT[®]

Demineralized Fibers



DEMINERALIZED FIBERS

Volume

1.0 cc

Order No. BL-1000-001

Freeze Dried

2.5 cc

Order No. BL-1000-002

Freeze Dried

5.0 cc

Order No. BL-1000-003

Freeze Dried

10 cc

Order No. BL-1000-004

Freeze Dried

1864 Concert Drive
Virginia Beach, VA 23453

www.LifeNetHealth.org

 For Use: Global

The LifeNet Health and ReadiGraft logos are registered trademark of LifeNet Health.
Allowash XG are registered trademarks of LifeNet Health.

©2015 LifeNet Health. All rights reserved.

68-60-024-01 .01

 **LifeNet Health[®]**
Saving Lives. Restoring Health. Giving Hope.

Phone: US 1-888-847-7831 • OUS 1-757-464-4761 x2000 • Fax: US 1-888-847-7832 or OUS 1-757-464-6439

www.LifeNetHealth.org