


FlexiGRAFT®

Demineralized Sponge




FlexiGRAFT® Demineralized Sponge is a demineralized cancellous matrix comprised of 100% cancellous bone that maintains natural bone architecture and acts as a conduit for vascular and cellular in-growth.

Brand	<ul style="list-style-type: none"> FlexiGRAFT® Demineralized Sponge
Clinical Applications	<ul style="list-style-type: none"> Any surgical application that requires bone void filler
Features And Benefits	<ul style="list-style-type: none"> Natural Osteoconductive Scaffold: The natural interconnected porous 3D architecture and increased surface area of demineralized cancellous provides an optimal scaffold for cellular attachment, and proliferation. Bioactive: Contains exposed natural growth factors with Osteoinductive potential. No Carrier: Contains no fillers or carriers commonly used in bone putties and gels. Hydrophilic: Naturally absorbs and retains bioactive fluids like blood, Platelet Rich Plasma (PRP), and Bone Marrow Aspirate (BMA). Biocompatible: FlexiGraft Demineralized Sponge is comprised of 100% natural bone matrix that is remodeled naturally during the bone healing process. Sterile: Sterilized using proprietary and patented Allowash® XG technology. This technology provides a sterility assurance level of 10⁻⁶, without compromising the grafts inherent osteoconductivity or osteoinductive potential Malleable, Elastic, Compressible: When re-hydrated, the product produces a sponge-like effect which allows the graft to be shaped and compressed to fit in and around a variety of bony defects.
Alternative Products	<ul style="list-style-type: none"> FlexiGraft Demineralized Fibers Cancellous Cortical Cancellous



FlexiGRAFT[®] Demineralized Sponge



CUBE

Size



8 x 8 x 8 mm
Order No. BL-1100-001
Freeze Dried



10 x 10 x 10 mm
Order No. BL-1100-003
Freeze Dried



12 x 12 x 12 mm
Order No. BL-1100-002
Freeze Dried



STRIP

Size



10 x 20 x 2 mm
Order No. BL-1300-001
Freeze Dried



15 x 40 x 2 mm
Order No. BL-1300-002
Freeze Dried



20 x 25 x 6 mm
Order No. BL-1300-003
Freeze Dried



10 x 20 x 8 mm
Order No. BL-1300-004
Freeze Dried



CHIPS

Grind Size Volume

1-4 mm 1.0 cc
Order No. BL-1200-001
Freeze Dried

1-4 mm 2.5 cc
Order No. BL-1200-002
Freeze Dried

1-4 mm 5.0 cc
Order No. BL-1200-003
Freeze Dried

1864 Concert Drive
Virginia Beach, VA 23453

www.AccessLifeNetHealth.org

The LifeNet Health and FlexiGraft logos are registered trademarks of LifeNet Health. ©2014 LifeNet Health, Virginia Beach, VA. All rights reserved.

68-60-026-01.00

