



FlexiGRAFT[®]

Demineralized Sponge

Clinical Overview

FlexiGraft Demineralized Sponge is a demineralized cancellous matrix comprised of 100% cancellous bone that maintains natural bone architecture and acts as a conduit for vascular and cellular in-growth.

Applications

- Any surgical application that requires bone void filler

Why Use

- **Compressible:** After re-hydration, graft can be compressed to fill tight spaces and will expand to fill the defect site after implantation
- **Bioactive:** Contains exposed natural growth factors with osteoinductive potential
- **Osteoconductive:** Natural bone matrix facilitates cell attachment and proliferation, and vascular in-growth
- **Sterile:** Sterilized using proprietary and patented Allowash XG[®] technology which provides a sterility assurance level of 10⁻⁶, without compromising the graft's inherent osteoconductive properties or osteoinductive potential¹
- **100% Human Bone:** Can remodel with the patient's own tissue during the healing process
- **Absorbent:** Absorbs and retains bioactive fluids like blood, platelet rich plasma (PRP), and bone marrow aspirate (BMA)
- **Versatile:** Available as strips, cubes, or chips in multiple sizes and volumes to meet surgical needs

References

1. Eisenlohr LM. "Allograft Tissue Sterilization Using Allowash XG[®]". 2007 Bio-Implants Brief.





Cube Freeze-dried, 5 year shelf life		Strip Freeze-dried, 5 year shelf life		Chips Freeze-dried, 3 year shelf life	
Size	Order Code	Size	Order Code	Size	Order Code
8 x 8 x 8 mm	BL-1100-001	10 x 20 x 2 mm	BL-1300-001	1-4 mm, 1.0 cc	BL-1200-001
10 x 10 x 10 mm	BL-1100-003	15 x 40 x 2 mm	BL-1300-002	1-4 mm, 2.5 cc	BL-1200-002
12 x 12 x 12 mm	BL-1100-002	20 x 25 x 6 mm	BL-1300-003	1-4 mm, 5.0 cc	BL-1200-003
		10 x 20 x 8 mm	BL-1300-004		

68-60-026-01.01

The LifeNet Health logo is a registered trademark of LifeNet Health.
FlexiGraft and Allowash XG are registered trademarks of LifeNet Health.
©2017 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu