

OraGRAFT®

Demineralized Large Particle



The OraGRAFT® product line is comprised of hard and soft allograft bio-implants generally used for sinus augmentation, craniofacial reconstruction, and correction of periodontal and ridge defects.

Brand

- OraGRAFT®

Clinical Applications

- Ridge Maintenance
- Extraction Site Preservation
- Sinus Grafting
- Craniofacial Reconstruction

Features And Benefits

- **Increased Safety** – Sterile with an SAL 10^{-6}
- **Promotes rapid healing** – Osteoinductive graft utilizing PAD technology that targets optimal residual calcium level 1-4%
- **Convenient Storage** – Ambient Temperature with 3 year shelf life (cortical), 5 year (cancellous)
- **Versatile** – Available as cortical or cancellous for variety of uses, larger particles hold space and allow for rapid vascular ingrowth



OraGRAFT®

Demineralized Large Particle



DEMINERALIZED LARGE PARTICLE

Demineralized Cortical 1 - 4 mm particles

Size	Freeze-Dried
5.0 cc	DGC5

Demineralized Cancellous 1 - 8 mm particles

Size	Freeze-Dried
5.0 cc	DCAN5

1864 Concert Drive
Virginia Beach, VA 23453
1-888-847-7831

The LifeNet Health and OraGraft logos are registered trademarks of LifeNet Health.
©2014 LifeNet Health, Virginia Beach, VA. All rights reserved.
68-60-074 .00

www.AccessLifeNetHealth.org

