

OrACELL[®]

Acellular Dermis



OrACELL[®] is an acellular dermal matrix used for soft tissue reconstruction, in areas such as periodontal defects and ridge preservation. It has $\geq 97\%$ DNA removal, comes ready-to-use right out of the package, can be stored at room temperature, and has an SAL (Sterility Assurance Level) of 10^{-6} , which is equal to medical device grade sterility.

Brand

- OrACELL[®]

Clinical Applications

- Periodontal Defects
- Ridge Preservation
- Soft Tissue Correction

Features And Benefits

- **Increased Safety** – Sterile with an SAL 10^{-6}
- **Convenient Storage** – Ambient Temperature
- **Biocompatible** – $\geq 97\%$ of DNA removed while retaining native growth factors, collagen and elastin
- **Ready to use** – No rehydration needed
- **Protects Graft Site** – Cell-occlusive natural membrane that protects grafted site from fluid and bacterial contamination
- **Effective Site Remodeling** – Graft resorbs through normal site remodeling in 4–6 months following placement



OrACELL[®]

Acellular Dermis



ACELLULAR DERMIS

Thickness = 0.76 - 1.25 mm

Size	Room Temperature
1.0 x 1.0 cm	OCELL 150
1.0 x 4.0 cm	OCELL 151
1.5 x 2.0 cm	OCELL100
2.0 x 4.0 cm	OCELL101

Thickness = 1.26 - 1.75 mm

Size	Room Temperature
1.0 x 1.0 cm	OCELL 250
1.0 x 4.0 cm	OCELL 251
1.5 x 2.0 cm	OCELL200
2.0 x 4.0 cm	OCELL201

1864 Concert Drive
Virginia Beach, VA 23453

The LifeNet Health logo is a registered trademark of LifeNet Health, Inc.
Oracell is a registered trademark of LifeNet Health, Inc., Virginia Beach, VA.

©2014 LifeNet Health. All rights reserved.

68-60-084-02 .00

www.AccessLifeNetHealth.org

