

# FlexiGRAFT®

## KinetiGRAFT®



FlexiGraft KinetiGraft is a pre-sutured tendon construct that is indicated for soft tissue approximation and/or ligation. It can be used for ACL and PCL Reconstruction procedures.

Brand	<ul style="list-style-type: none"> <li>FlexiGraft KinetiGraft</li> </ul>
Clinical Applications	<ul style="list-style-type: none"> <li>ACL and PCL Reconstruction</li> </ul>
Features And Benefits	<ul style="list-style-type: none"> <li><b>Convenience:</b> No prep time or timely harvest of autograft required, no training for staff necessary. An out-of-the-box option for OR efficiency.</li> <li><b>Consistency:</b> Trained technicians suturing graft for construct consistency.</li> <li><b>Compatible:</b> Construct was designed by Arthrex® to be used with their specialized tools, fixation and guides that have already been developed and are in use by surgeons today.</li> <li><b>Patient-friendly:</b> Construct eliminates donor site morbidity and associated pain from autograft harvest. This makes the procedure less invasive and potentially decreases OR time. Less OR time can mean less time under anesthesia and less tourniquet time.</li> <li><b>Sterile:</b> KinetiGraft is sterilized using proprietary and patented Allowash XG® technology. This provides a sterility assurance level of 10<sup>-6</sup>, without compromising the construct's inherent biomechanical properties.</li> </ul>
Alternative Products	<ul style="list-style-type: none"> <li>Tibialis &amp; Hamstring Tendons</li> <li>Peroneus Longus Tendons</li> <li>Patellar &amp; Achilles Tendons</li> <li>GraftLink® Construct</li> </ul>

# FlexiGRAFT®

## KinetiGRAFT®



### KinetiGRAFT

Diameter 9.5 - 10.5 mm

Length: 80 - 90 mm

Order No. FKG10

Frozen

1864 Concert Drive  
Virginia Beach, VA 23453

[www.LifeNetHealth.org](http://www.LifeNetHealth.org)

**US** For Use: US only

The LifeNet Health logo, FlexiGraft, KinetiGraft and Allowash XG are registered trademarks of LifeNet Health, Inc., Virginia Beach, VA. GraftLink is a registered trademark of Arthrex, Inc.

©2015 LifeNet Health. All rights reserved.

68-60-III.00

 **LifeNet Health®**  
Saving Lives. Restoring Health. Giving Hope.

Phone: US 1-888-847-7831 • OUS 1-757-464-4761 x2000 • Fax: US 1-888-847-7832 or OUS 1-757-464-6439

[www.LifeNetHealth.org](http://www.LifeNetHealth.org)