



# FlexiGRAFT<sup>®</sup>

## Pre-Sutured Lateral Ankle Tendon

### Clinical Overview

FlexiGraft Pre-Sutured Lateral Ankle Tendon Constructs offer out-of-the-box convenience, eliminate donor site morbidity and pain issues due to autograft harvests and are medical device grade sterile (SAL 10<sup>-6</sup>)

### Applications

- Lateral Ankle Reconstruction

### Why Use

- **Convenience:** No prep time or timely harvest of autograft required, no training for staff necessary. An out-of-the-box option for OR efficiency
- **Pre-Sized:** Construct designed to a diameter of 4 to 5 mm and length of 150 to 160 mm – eliminating the need for trimming by the surgeon's staff
- **Consistency:** Trained technicians suturing graft for construct consistency. Removes the variability between surgical assistants or physician assistants preparing the tendon
- **Patient-Friendly:** Construct eliminates donor site morbidity and associated pain from autograft harvest. This makes the procedure less invasive and likely decreases OR time. Less OR time can mean less time under anesthesia and less tourniquet time
- **Sterile:** Pre-Sutured Lateral Ankle Tendon is sterilized using proprietary and patented Allowash XG<sup>®</sup> technology. This provides a sterility assurance level (SAL) of 10<sup>-6</sup>, without compromising the construct's inherent biomechanical properties





## Pre-Sutured Lateral Ankle Tendon

Frozen Storage/3 year shelf life

Diameters	Length	Order Code
4.0, 4.5 and 5.0 mm	150 - 160 mm	FPSST

### How do we size Pre-Sutured Lateral Ankle Construct?

Length is measured to be 150 - 160 mm without tension.

Diameter is measured single strand by pulling the construct through a tendon sizer with modest pressure. The recorded diameter is the smallest channel the tendon will pass through.

*\*This graft may not be available in all countries.*

68-60-112-01.02

The LifeNet Health logo is a registered trademark of LifeNet Health.  
FlexiGraft and Allowash XG are registered trademarks of LifeNet Health.  
©2019 LifeNet Health, Virginia Beach, VA. All rights reserved.



**North America**  
1.888.847.7831  
orders@lifenethealth.org

**Europe**  
+ 43 1 375002710  
eu\_orders@lifenethealth.eu

**Latin America** □ **Asia** □ **Middle East**  
1.757.464.4761 ext. 2000  
internat.orders@lifenethealth.org

**LifeNetHealth.org**  
**LifeNetHealth.eu**