


# Costal Cartilage



Costal Cartilage is sterile ( $10^{-6}$ ) hyaline cartilage from the rib. Cleaned and sterilized using the patented Allowash XG<sup>®</sup> technology, this grafts biomechanical and biochemical properties are maintained, making it ideal for nasal and ear reconstruction.

Product Name	<ul style="list-style-type: none"><li>• Costal Cartilage</li></ul>
Clinical Applications	<ul style="list-style-type: none"><li>• Nasal Reconstruction</li><li>• Ear Reconstruction</li></ul>
Features And Benefits	<ul style="list-style-type: none"><li>• <b>Sterile:</b> Sterilized using proprietary and patented Allowash XG technology. This technology provides a sterility assurance level of <math>10^{-6}</math>, without compromising the grafts inherent properties</li><li>• <b>Biocompatible:</b> Costal Cartilage is comprised of 100% hyaline cartilage that is remodeled naturally during the healing process</li><li>• <b>Supportive:</b> Provides biomechanical strength requirements for intended use</li><li>• <b>Convenient:</b> Allows physician to select appropriate graft for clinical procedure</li><li>• <b>Malleable:</b> Easy to trim and shape to perfectly fit surgical need</li></ul>

# Costal Cartilage



Costal Cartilage

Size	Preservation Method	Order Number
L: 30 - 59 mm; D: $\geq$ 7 mm; W: $\geq$ 7 mm	Frozen	CCART30
L: 60 - 90 mm; D: $\geq$ 7 mm; W: $\geq$ 7 mm	Frozen	CCART60

1864 Concert Drive  
Virginia Beach, VA 23453

[www.LifeNetHealth.org](http://www.LifeNetHealth.org)

The LifeNet Health logo and Allowash XG are registered trademarks of LifeNet Health.

©2016 LifeNet Health. All rights reserved.

68-60-113 .01

