



# VertiGRAFT® Spacers

---

## Clinical Overview

- Solid one piece cortical spacer designed to provide biomechanical strength for spinal fusion inter-body and discectomy procedures

---

## Applications


- Spinal fusion inter-body and discectomy procedures

---

## Why Use

- **Strength:** Cortical Construction
- **Stability:** Textured surface area to reduce migration
- **Variety:** Multiple footprints in both parallel and lordotic configurations to better match patient anatomy
- **Sterile:** Sterilized using proprietary and patented Allowash XG® technology. This technology provides a sterility assurance level of  $10^{-6}$ , without compromising the implant's inherent osteoconductive properties<sup>1</sup>
- **Convenience:** Freeze-dried for room temperature storage

1. Eisenlohr LM. "Allograft Tissue Sterilization Using Allowash XG®" 2007 Bio-Implants Brief.



## Fibular Lordic

Shape	Freeze Dried	Diameter	Anterior Height H(A)*	Posterior Height H(P)*
Lordotic	VFIBW5	8-15 mm	8 mm	5 mm
Lordotic	VFIBW7	8-15 mm	10 mm	7 mm
Lordotic	VFIBW9	8-15 mm	12 mm	9 mm
Lordotic	VFIBW11	8-15 mm	14 mm	11 mm

## Fibular Parallel

Shape	Freeze Dried	Diameter	Anterior Height H(A)*
Parallel	VFIBCS16	8-15 mm	7 mm
Parallel	VFIBCS17	8-15 mm	8 mm
Parallel	VFIBCS18	8-15 mm	9 mm
Parallel	VFIBCS19	8-15 mm	10 mm
Parallel	VFIBCS110	8-15 mm	11 mm
Parallel	VFIBCS112	8-15 mm	12 mm
Parallel	VFIBCS114	8-15 mm	13 mm

## Humeral Lordic

Shape	Freeze Dried	Diameter	Anterior Height H(A)*	Posterior Height H(P)*
Lordotic	VHUMW11	18-22 mm	14 mm	11 mm
Lordotic	VHUMW13	18-22 mm	16 mm	13 mm
Lordotic	VHUMW15	18-22 mm	18 mm	15 mm
Lordotic	VHUMW17	18-22 mm	20 mm	17 mm

## Humeral Parallel

Shape	Freeze Dried	Diameter	Anterior Height H(A)*
Parallel	VHUMCS8	18-22 mm	8 mm
Parallel	VHUMCS10	18-22 mm	10 mm
Parallel	VHUMCS12	18-22 mm	12 mm
Parallel	VHUMCS14	18-22 mm	14 mm
Parallel	VHUMCS16	18-22 mm	16 mm

\*Nominal Measurements

LifeNet Health ■ 1864 Concert Drive ■ Virginia Beach, VA 23453

68-60-128.00 REV.00

The LifeNet Health logo is a registered trademark of LifeNet Health, Inc.  
 Allowash XG and VertiGraft logos are registered trademarks of LifeNet Health, Inc.  
 ©2016 LifeNet Health. All rights reserved.

