



FlexiGRAFT[®]

Achilles Tendons

Clinical Overview

FlexiGraft Achilles Tendons can be used for ACL and PCL Reconstruction procedures


Applications

- ACL and PCL Reconstruction
- Multiligament Repairs
- Achilles Tendon Repair
- Other Soft Tissue Repairs

Why Use

- **Convenience:** No timely harvest required. Pre-shaped grafts eliminate prep time and further increase operating room efficiency
- **Patient-Friendly:** Eliminates donor site morbidity and associated pain from autograft recovery. This makes the procedure less invasive and potentially decreases OR time. Less OR time can mean less time under anesthesia and less tourniquet time
- **Sterile:** Achilles tendons are sterilized using proprietary and patented Allowash XG[®] technology. This provides a sterility assurance level (SAL) of 10⁻⁶, without compromising the construct's inherent biomechanical properties





Achilles Tendons

Frozen Storage/5 year shelf life

Description	Tendon Length (TL)	Bone Block Size	Order Code
Without Bone Block	160 - 260 mm	n/a	FAT
With Bone Block	160 - 260 mm	width = 10 - 26 mm (W) length = 25 - 30 mm (L) thickness = 10 mm (T)	FATB

Description	Size	Order Code
Pre-Shaped Bone Block	160 - 260 mm 10 mm diameter dowel length 25 - 30 mm	FATB10
	160 - 260 mm 11 mm diameter dowel length 25 - 30 mm	FATB11

68-60-108.01

The LifeNet Health logo is a registered trademark of LifeNet Health.
FlexiGraft and Allowash XG are registered trademarks of LifeNet Health.
©2017 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu