



ArthroFLEX[®]

Acellular Dermal Matrix

Clinical Overview

- ArthroFlex is an acellular dermal matrix used to supplement soft tissue repairs such as the Achilles tendon and rotator cuff. It has $\geq 97\%$ DNA removal, comes ready-to-use right out of the package, can be stored at room temperature, and has a SAL (Sterility Assurance Level) of 10^{-6} , which is equal to medical device grade sterility

Applications

- Augment repairs of the Rotator Cuff
- Augment repairs of the Achilles, Biceps, Quad and Patellar tendons
- Reinforce peroneal/tibials tendons and flexor/extensor tendons of hand and foot
- Reinforce capsular repairs
- Suture reinforcement of soft tissue repairs

Why Use

- **Increased Safety:** Sterile with a SAL of 10^{-6} means lower risk of disease transmission
- **Decellularized:** Uses Matracell[®] Technology to remove $\geq 97\%$ DNA and cellular remnants to decrease likelihood of immune response¹
- **Promotes Rapid Healing:** Biocompatible and retains natural growth factors, collagen and elastin¹
- **Convenient:** Preservon[®] technology allows graft to be stored fully hydrated at room temperature so it can be used right out of the package, without lengthy OR reconstitution processes
- **Versatility:** Available in multiple sizes and thicknesses

References

1. Data on file at LifeNet Health

Arthrex[®]

LifeNet Health[®]
Saving Lives. Restoring Health. Giving Hope.

North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu



ARTHOFLEX ACELLULAR DERMAL MATRIX

Room Temperature Storage / 3 Year Shelf Life

Size	Thickness	Order Code
35 x 35 mm	1.5 mm (1.26 - 1.74)	AFLEX100
40 x 70 mm	1.5 mm (1.26 - 1.74)	AFLEX101
50 x 90 mm	1.5 mm (1.26 - 1.74)	AFLEX103
15 x 140 mm	1.5 mm (1.26 - 1.74)	AFLEX150
35 x 35 mm	2.0 mm (1.26 - 1.74)	AFLEX200
40 x 70 mm	2.0 mm (1.26 - 1.74)	AFLEX201
40 x 70 mm	3.0 mm (2.5 - 3.5)	AFLEX301
40 x 70 mm	1.0 mm (1.26 - 1.74)	AFLEX400
40 x 70 mm	1.0 mm (1.26 - 1.74)	AFLEX401
30 x 40 mm	0.5 mm (0.3 - 1.0)	AFLEX500
BioWasher™ 10 x 14 mm, 2 pack	2.0 mm (1.76 - 2.24)	AFLEX822*

* AFLEX822 has an 18 Month Shelf Life

68-60-081-02.00

The LifeNet Health logo, ArthroFlex, Preservon and Matracell are registered trademarks of LifeNet Health. BioWasher is a trademark of LifeNet Health. ©2017 LifeNet Health. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu