



# ArthroFLEX

## BioWasher™

### Acellular Dermal Matrix

#### Clinical Overview

- ArthroFlex BioWasher is a smaller sized acellular dermal matrix intended for suture reinforcement in soft tissue repairs, such as the rotator cuff


#### Applications

- Suture reinforcement of soft tissue repairs
- Arthroscopic Rotator Cuff repair
- Plantar Plate repair

#### Why Use

- **Smaller Size:** Easier to implant in arthroscopic procedures, easily passes through surgical portal
- **Increased Safety:** Sterile with a SAL of  $10^{-6}$  means lower risk of disease transmission
- **Decellularized:** Uses Matracell® Technology to remove DNA and cellular remnants to decrease likelihood of immune response<sup>1</sup>





## BioWasher

Room Temperature Storage / 18 Month Shelf Life

Size	Thickness	Order Code
10 x 14 mm rounded rectangle, 2 per package	2.0 mm (1.76 - 2.24)	AFLEX822

\*

68-60-093-01.02

The LifeNet Health logo is a registered trademark of LifeNet Health, Inc. Arthroflex and Matracell are registered trademarks of LifeNet Health, Inc. ©2017 LifeNet Health. All rights reserved.



**North America**  
1.888.847.7831  
orders@lifenethealth.org

**Europe**  
+ 43 1 375002710  
eu\_orders@lifenethealth.eu

**Latin America** □ **Asia** □ **Middle East**  
1.757.464.4761 ext. 2000  
internat.orders@lifenethealth.org

[www.LifeNetHealth.org](http://www.LifeNetHealth.org)  
[www.LifeNetHealth.eu](http://www.LifeNetHealth.eu)