



OraGRAFT®

Cancellous Blocks and Cubes

Clinical Overview

The OraGraft product line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects

Applications

- Ridge Maintenance (Horizontal & Vertical Augmentation)
 - Craniofacial Reconstruction
-

Why Use

- **Increased Safety:** Sterile with an SAL 10^{-6}
- **Convenient Storage:** Ambient Temperature with 5 year shelf life
- **Ease of Use:** Can be affixed to surgical site using standard bone screws
- **Rapid Remodeling:** Solid cancellous bone allows for space maintenance and rapid remodeling due to trabecular architecture
- **Match to Defect:** Easy to trim, can be shaped using standard instrumentation





Cancellous Blocks and Cubes

Ambient Temperature*/5 year shelf life

Size	Freeze-Dried
10 x 10 x 10 mm	CANCUBE1
15 x 15 x 15 mm	CCUBE-01
15 x 15 x 8 mm	CCUBE-02
15 x 30 x 8 mm	CCUBE-03

*While ambient temperature has not been defined by AATB, LifeNet Health maintains its lyophilized tissue at 10-37°C

68-60-079 .02

The LifeNet Health logo is a registered trademark of LifeNet Health.
OraGraft is a registered trademark of LifeNet Health.
©2017 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu