



OraGRAFT[®]

Cancellous Mineralized Powder

Clinical Overview

The OraGraft product line is comprised of hard and soft allograft bio-implants generally used for sinus augmentation, craniofacial reconstruction, and correction of periodontal and ridge defects

Applications

- Periodontal Defects
- Ridge Maintenance
- Extraction Site Preservation
- Sinus Grafting
- Craniofacial Reconstruction

Why Use

- **Increased Safety:** Sterile with an SAL 10^{-6}
- **Convenient Storage:** Ambient Temperature with 3 year shelf life
- **Graft Incorporation:** Osteoconductive graft with an open trabecular structure allows for healing through “creeping substitution”





Cancellous Mineralized Powder

Ambient Temperature*/3 year shelf life

250 - 1000 microns

Size	Freeze-Dried
0.50 cc	OCAN-0.5A
1.00 cc	OCAN-1.0A
2.00 cc	OCAN-2.0A

1000 - 2000 microns

0.50 cc	OCAN-0.5B
1.00 cc	OCAN-1.0B
2.00 cc	OCAN-2.0B

*While ambient temperature has not been defined by AATB, LifeNet Health maintains its lyophilized tissue at 10-37°C

68-60-077 .01

The LifeNet Health logo is a registered trademark of LifeNet Health.
Oragraft is a registered trademark of LifeNet Health.
©2017 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu