



# FlexiGRAFT®

## Cannulated Revision Dowels

### Clinical Overview

FlexiGraft Revision Dowels are 100% cancellous bone with a thin cortical plate for use in backfilling tunnels created during ACL/PCL procedures

### Applications

- Used to fill tunnels created for revision ACL/PCL procedure
- Other procedures/bone voids that require a dowel to backfill autologous sites

### Why Use

- **Biocompatible:** Comprised of 100% human cancellous bone with a thin cortical plate that is remodeled naturally during the bone healing process
- **Sterile:** Sterilized using proprietary and patented Allowash XG® technology. This technology provides a sterility assurance level (SAL) of  $10^{-6}$ , without compromising the implant's inherent osteoconductive properties
- **Fully Hydrated, Ambient Storage Available:** This implant features Preservon®, a proprietary, glycerol-based preservation technology that allows allograft bio-implants to be stored in a fully hydrated state at ambient temperature. Preservon eliminates lengthy rehydrating times, and does not require freezer storage
- **Convenience:** Fully cannulated with chamfered tip for easy insertion/implantation

### Reference

1. Hunt, D.G. & Chandler, C & Ulrich, D.A. & Poska, R & Montero, A. (2013). USP council of experts USP controlled room temperature range expansion. Pharmacopeial Forum. 39.





Cannulation Diameter = 2.49 mm

Standard Length Dowels 25-29 mm	
Ambient Storage*/5 year shelf life	
Order Code	Sizing
PCD9	9 mm diameter
PCD10	10 mm diameter
PCD11	11 mm diameter
PCD12	12 mm diameter
PCD13	13 mm diameter
PCD14	14 mm diameter
PCD16	16 mm diameter
PCD18	18 mm diameter

Extra Long Dowels 30-35 mm	
Ambient Storage*/5 year shelf life	
Order Code	Sizing
PCDXL9	9 mm diameter
PCDXL10	10 mm diameter
PCDXL11	11 mm diameter
PCDXL12	12 mm diameter
PCDXL13	13 mm diameter
PCDXL14	14 mm diameter
PCDXL16	16 mm diameter
PCDXL18	18 mm diameter

\*While ambient room temperature has not been defined by the AATB, the expanded controlled room temperature range suggested by the U.S. Pharmacopeial (USP) stimulus document is 2°C to 30°C with excursions up to 40°C.<sup>1</sup>

68-60-101-01 .03

LifeNet Health, the LifeNet Health logo, FlexiGraft, Allowash XG and Preservon are registered trademarks of LifeNet Health.  
©2019 LifeNet Health, Virginia Beach, VA. All rights reserved.



**North America**  
1.888.847.7831  
orders@lifenethealth.org

**Europe**  
+ 43 1 375002710  
eu\_orders@lifenethealth.eu

**Latin America** □ **Asia** □ **Middle East**  
1.757.464.4761 ext. 2000  
internat.orders@lifenethealth.org

[www.LifeNetHealth.org](http://www.LifeNetHealth.org)  
[www.LifeNetHealth.eu](http://www.LifeNetHealth.eu)