



# OraGRAFT<sup>®</sup>

## Cortical Demineralized Powder

### Clinical Overview

The OraGraft product line is comprised of hard and soft allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects

### Applications

- Periodontal Defects
- Ridge Maintenance
- Extraction Site Preservation
- Sinus Grafting
- Craniofacial Reconstruction

### Why Use

- **Highest Safety Level:** Sterile with an SAL of  $10^{-6}$
- **Convenient Storage:** Ambient temperature with a 3 year shelf life
- **Promotes Rapid Healing:** Osteoinductive graft utilizing PAD<sup>®</sup> technology that targets optimal residual calcium level 1-4%
- **Versatile:** Available in 5 size configurations





## Cortical Demineralized Powder

Ambient Temperature\*/3 year shelf life

250 - 1000 microns

Size	Freeze-Dried
0.25 cc	DGC1/20
0.50 cc	DGC1/10
0.70 cc	DGC1/8
1.20 cc	DGC1/4
2.50 cc	DGC

\*While ambient temperature has not been defined by AATB, LifeNet Health maintains its lyophilized tissue at 10-37°C

68-60-078 .02

The LifeNet Health logo is a registered trademark of LifeNet Health.  
Oragraft is a registered trademark of LifeNet Health.  
©2017 LifeNet Health, Virginia Beach, VA. All rights reserved.



**North America**  
1.888.847.7831  
orders@lifenethealth.org

**Europe**  
+ 43 1 375002710  
eu\_orders@lifenethealth.eu

**Latin America** □ **Asia** □ **Middle East**  
1.757.464.4761 ext. 2000  
internat.orders@lifenethealth.org

**www.LifeNetHealth.org**  
**www.LifeNetHealth.eu**