

ReadiGRAFT®

Demineralized Chips and Particulate



ReadiGRAFT®
Demineralized Cortical Particulate and partially demineralized cancellous and cortical cancellous chips.

Brand

- ReadiGRAFT®

Clinical Applications

- Any surgical application that requires bone void filler

Features And Benefits

- **Natural Osteoconductive Scaffold:** provide an optimal scaffold for cellular attachment, proliferation, and vascular in-growth.
- **Hydrophilic:** Absorbs and retains bioactive fluids like blood, Platelet Rich Plasma (PRP), and Bone Marrow Aspirate (BMA).
- **Sterile:** The implant is sterilized using proprietary and patented Allowash XG® technology. This technology provides a sterility assurance level of 10^{-6} , without compromising the implant's inherent osteoconductivity or osteoinductive potential.
- **Osteoinductive Potential:** Demineralized using a proprietary demineralization process that preserves naturally occurring growth factors.

Alternative Products

- Cancellous Chips
- Cortical Cancellous
- Femoral Head (Grinder)



ReadiGRAFT®

Demineralized Chips and Particulate



CANCELLOUS CHIPS

Grind Size	Volume	Freeze Dried Order No.
1-8 mm	40 cc	DCAN40

CORTICAL CHIPS

Grind Size	Volume	Freeze Dried Order No.
125-1000 microns	20 cc	DGC20
125-1000 microns	40 cc	DGC40

CORTICAL CANCELLOUS CHIPS

Grind Size	Volume	Freeze Dried Order No.
1-8 mm	20 cc	DCC1/4

1864 Concert Drive
Virginia Beach, VA 23453
757-464-4761 x2000
Fax: 757-464-6439

www.AccessLifeNetHealth.org

The LifeNet Health and Readigraft logos are registered trademarks of LifeNet Health.
©2012 LifeNet Health, Virginia Beach, VA. All rights reserved.
68-50-252.01

