



# MatriGRAFT<sup>®</sup>

## Fibular Wedge

### Clinical Overview

Natural shaped, parallel fibular wedge, designed to provide immediate structural support and restore segmental bone loss

### Applications

- Anterior Cervical Fusion
- Trauma

### Why Use

- **Osteoconductive:** Natural bone matrix facilitates cell attachment and proliferation, and vascular in-growth
- **Structural:** All-cortical design provides immediate structural support
- **100% Human Bone:** Will remodel alongside patient's own tissue during the healing process
- **Sterile:** Sterilized using patented and proprietary Allowash XG<sup>®</sup> technology which provides a sterility assurance level (SAL) of 10<sup>-6</sup>, without compromising the graft's inherent osteoconductive properties<sup>1</sup>
- **Convenient:** Implant is pre-sized to fit a variety of applications and minimize prep time in the operating room

### References

1. Eisenlohr LM. "Allograft Tissue Sterilization Using Allowash XG." 2007 Bio-Implants Brief.





## MatriGraft Fibular Wedge

Diameter: 8 - 15 mm ▪ Cortical Wall Thickness: 2 mm or greater

Freeze-Dried Code	Frozen Code	Height
FIBW0.6		6 mm
FIBW0.8		8 mm
FIBW1		10 mm
FIBW2		12 mm
FIBW4		14 mm
FIBW6		16 mm
FIBW8	FFIBW8	18 mm

68-60-037-01.00

The LifeNet Health logo is a registered trademark of LifeNet Health.  
 MatriGraft and Allowash XG are registered trademarks of LifeNet Health.  
 ©2017 LifeNet Health, Virginia Beach, VA. All rights reserved.



**North America**  
 1.888.847.7831  
[orders@lifenethealth.org](mailto:orders@lifenethealth.org)

**Europe**  
 + 43 1 375002710  
[eu\\_orders@lifenethealth.eu](mailto:eu_orders@lifenethealth.eu)

**Latin America** □ **Asia** □ **Middle East**  
 1.757.464.4761 ext. 2000  
[internat.orders@lifenethealth.org](mailto:internat.orders@lifenethealth.org)

[www.LifeNetHealth.org](http://www.LifeNetHealth.org)  
[www.LifeNetHealth.eu](http://www.LifeNetHealth.eu)