



FlexiGRAFT[®]

Patellar Tendons

Clinical Overview

FlexiGraft patellar tendons can be used for ACL or PCL Reconstruction procedures

Applications

- ACL Reconstruction
- PCL Reconstruction


Why Use

- **Convenience:** No timely harvest required. Pre-shaped grafts eliminate prep time and further increase operating room efficiency
- **Appropriate for ACL:** Considered Gold Standard
- **Patient-Friendly:** Eliminates donor site morbidity and associated pain from autograft recovery. This makes the procedure less invasive and potentially decreases OR time. Less OR time can mean less time under anesthesia and less tourniquet time¹
- **Sterile:** Patellar tendons are sterilized using proprietary and patented Allowash XG[®] technology. This provides a sterility assurance level (SAL) of 10⁻⁶, without compromising the construct's inherent biomechanical properties

References

1. Oro et al. Autograft Versus Allograft: An Economic Cost Comparison of Anterior Cruciate Ligament Reconstruction. Arthroscopy 2011; 27(9):1219-1225.





Patellar Tendons

Frozen Storage/5 year shelf life

Description	Tendon Length (TL) Tendon length is measured as length between bone blocks	Size info	Order Code
Bisected Patellar Tendon	30 - 55 mm	Tendon insertion width (TI) = 10 - 25 mm Bone block ≥ 10 x 12 x 25 mm (W x T x L)	FBPL
Bisected Patellar Tendon with Small Block	30 - 55 mm	Tendon insertion width (TI) = 10 - 25 mm Bone block ≥ 10 x 10 x 20 mm (W x T x L)	FBPLSB
Whole Patellar Tendon	30 - 55 mm	Tendon insertion width (TI) = 20 - 35 mm Bone block ≥ 20 x 12 x 25 mm (W x T x L)	FWPL
Whole Patellar Tendon with Quad	30 - 55 mm	Tendon insertion width (TI) = 20 - 35 mm Bone block ≥ 20 x 12 x 25 mm (W x T x L) Quad tendon > 60 mm	FWPLQ
Whole Patellar Tendon with Small Block	30 - 55 mm	Tendon insertion width (TI) = 20 - 35 mm Bone block ≥ 20 x 10 x 20 mm (W x T x L)	FWPLSB
Pre-Shaped Patellar Tendon	35 - 50 mm	10 mm diameter dowel Tibial Dowel length 25 - 30 mm	FPL10
Pre-Shaped Patellar Tendon	35 - 50 mm	11 mm diameter dowel Tibial Dowel length = 25 - 30 mm	FPL11

68-60-109 .02

The LifeNet Health logo is a registered trademark of LifeNet Health.
FlexiGraft and Allowash XG are registered trademarks of LifeNet Health.
©2018 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu