

FlexiGRAFT®

Pre-Sutured Lateral Ankle Tendon



FlexiGRAFT Pre-Sutured Lateral Ankle Tendon Constructs offer out-of-the-box convenience, are compatible with the Lateral Ankle Implant Kit from Arthrex, eliminate donor site morbidity and pain issues due to autograft harvests and are medical device grade sterile (SAL 10⁻⁶).

Brand

- FlexiGraft

Clinical Applications

- Lateral Ankle Reconstruction

Features And Benefits

- **Convenience:** No prep time or timely harvest of autograft required, no training for staff necessary. An out-of-the-box option for OR efficiency.
- **Compatible:** Construct was designed by Arthrex to be used with their Lateral Ankle Implant Kit (AR-1675BC-CP). Specialized tools, fixation, and guides have already been developed and are in use by surgeons today.
- **Pre-sized:** Construct designed to a diameter of 4 to 5mm and length of 150 to 160mm - eliminating the need for trimming by the surgeon's staff.
- **Consistency:** Trained technicians suturing graft for construct consistency. Removes the variability between surgical assistants or physician assistants preparing the tendon.
- **Patient-Friendly:** Construct eliminates donor site morbidity and associated pain from autograft harvest. This makes the procedure less invasive and likely decreases OR time. Less OR time can mean less time under anesthesia and less tourniquet time.
- **Sterile:** Pre-Sutured Lateral Ankle Tendon is sterilized using proprietary and patented Allowash XG technology. This provides a sterility assurance level of 10⁻⁶, without compromising the construct's inherent biomechanical properties.

Alternative Products

- FSTP
- FST
- FGRACILIS



FlexiGRAFT[®]

Pre-Sutured Lateral Ankle Tendon



Pre-Sutured Lateral Ankle Tendon

Diameter: 4 to 5mm
Length: 155mm +/- 5mm
Order Code: FPSST
Frozen

1864 Concert Drive
Virginia Beach, VA 23453
Phone: 757-464-4761 x2000
Fax: 757-464-6439

www.AccessLifeNetHealth.org

The LifeNet Health and FlexiGraft logos are registered trademarks of LifeNet Health.
©2013 LifeNet Health, Virginia Beach, VA. All rights reserved.
68-60-048-01 Rev. 00

